## SPIS TREŚCI

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**INFRARED THERMOMETER**

**MODEL: TMB-COMPACT**

## PRODUCT INTRODUCTION, CLASSIFICATION AND COMPOSITIONE

This product measures the body temperature through gathering the infrared energy radiated from the forehead. It is easy, hygeian and accurate. The user only need to point the infrared probe to the forehead and press measure button, the body temperature will be shown on the LCD screen in 0.5 seconds. The product is widely used in home healthcare and medical institutes to measure body temperature. It is class B medical device with internal power supply and belongs to B-type application part. Its protection rate is IPX0. In EU classification, it is class IIa. It can not be used in any flammable gas and air, oxygen, or a combination of nitrous oxide gas. This product is mainly composed of infrared sensor, signal processer, button, buzzer, LCD screen and battery. To ensure you understand all safety information of the functions, please read carefully this manual before usage.

## BASIC OPERATION PRINCIPLE

Any object emits infrared radiation energy whose energy magnitude and wavelength is determined by its surface temperature. The wavelength from human body is 5~15  $\mu\text{m}$ . This product adopts the imported high precise infrared sensor to measure the wavelength released by human body to accurately measure the body temperature with precise calculation and various compensation correction.

## PRODUCT FEATURES

- With stable and reliable imported high precision infrared sensors
- Good performance and adaptability in the complex environment
- Patent probe design to ensure a more precise measurement


- One button operation for Body temperature and Object temperature
- Fever alarm, can set the alarm value
- Automatically memorize 32 measurement data.
- LCD screen, single color backlight
- Celsius(°C) and Fahrenheit (°F)for options
- Automatically power off to save

## TECHNICAL PARAMETERS

Temperature unit	°C/°F
Measurement method	Non-Contact
Effective distance	1cm-5cm
Measurement Mode	Body Mode: 32~42.9°C (89.6~109.2°F)
	Surface Mode: 0~100°C (32~212°F)
Measurement accuracy	Body Mode: Within (35-42)°C, ±0.2°C/0.4°F Beyond (35-42)°C, ±0.3°C/0.54°F
	Surface Mode ±1.0°C/°F
Display resolution	0.1°C/°F
Operating environment	15~35°C (59~95°F) RH<85%(moisture condensation)
Storage condition	-25~55°C (-13~-131°F) Rh < 90%(moisture condensation)
Power supply	DC 3V (2AAA 7#)
Power consumption	When off < 3 uA When measuring< 26 mA
Auto power-off	In 30 seconds
Unit Dimensions	43mm * 34 mm *155 mm(length x width x height)
Weight	76 g (without battery)

Memory	Save 32 sets measurement data
Electric quantity prompt	Low battery hint
Useful life	5 years

## SHAPE STRUCTURE


1. LCD display
2. Body temperature
3. Object temperature
4. Function button
5. Measure button
6. Infrared sensor
7. Label
8. Sensor cover
9. Handle
10. Battery cover

## 4.1 Display icon definition

	Show:battery is too low to work, please replace the battery Not show:enough battery
	Buzzer is off, no prompt sound
	body temperature measurement mode
	object temperature measurement mode
°C °F	Celsius degree Fahrenheit degree
	Temperature value
	Memorized measure data

## 4.2 The function button


Function button	Functional description
Setting menu button	Press to view the memorized measurement data, long press menu button to switch mode
 („+“)	Adjust the Settings of parameters,view measure, data upwards
 („-“)	Adjust the Settings of parameters,view measure, data downwards

## SETTINGS

This product has 4 function settings which is temperature units, sound indication switch, temperature alarm point, temperature offset and measure mode. Measure mode setting is set by the mode switch button. Other settings are set in the menu. The Settings menu table is as follows:

Menu	Function	“-“	+“	default value	Re- marks
F-1	Temperature alarm point	°C	°F	°C	
F-2	Temperature alarm point	Reduce to 0.1°C	Increase of 0.1°C	38°C	Scope of 37°C~42.9°C
F-3	Temperature offset	The downward migration of 0.1 °C	The upward migration of 0.1 °C	0.0°C	The range of -2° C ~ 2 °C
F-4	Indication switch	Off	On	On	

1. Battery
2. Body/Object temperature mode
3. Degree centigrade
4. Numerical display
5. Sound
6. Memory value data
7. Memory value symbol

## 1. Temperature measure mode setting

Turn on the device, the screen shows all the icons (Figure 1), start measuring temperature, press  $\Omega$  („+”) or  $\square$  („-”) to set the mode according to your requirement (Figure2 or 3).

Figure 2

Figure 3

Figure 4

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Figure 6

### 3. Temperature alarm point: F-2

Press twice the function button, F-2 shows on the screen, initial alarm point is 38.0° C (Figure 6), press  $\curvearrowright$  („+”) or  $\curvearrowleft$  („-”) to increase (or decrease) 0.1° C (Figure 7)

Figure 7

Figure 8

### 4. Temperature offset setting: F-3

Press thrice the function button, F-3 shows on the screen, initial value is 0.0°C, press  $\curvearrowright$  („+”) to shift upwards 0.1°C (Figure 8), press  $\curvearrowleft$  („-”) to shift downwards 0.0°C (Figure 9).

Figure 9

Figure 10

### 5. Sound indication switch setting: F-4

Press four times the function button, F-4 shows on the screen, press  $\curvearrowright$  („+”) or  $\curvearrowleft$  („-”) to turn on/off the sound indication (Figure 10, Figure 11), then press the function button again for confirm and exit.


Figure 11

Figure 12

## MEASUREMENT METHOD

### 1 Body temperature measurement.

- Press measure button to boot, all icons shows on the screen (Figure 12), single color as backlight, after beep it display the current measuring value (Figure 13). It is in the state of temperature measurement.
- Make sure the measure mode is body temperature mode (Figure 14).
- Aim the thermometer probe vertically to the forehead at the distance, 1-5 cm (Figure 15). Press measure button, after beep in 0.5s the measurement result shows on the screen.

Prompt  $32^{\circ}\text{C} \leq \text{measured body temperature} \leq 37.4^{\circ}\text{C}$ : normal  
 $37.5^{\circ}\text{C} \leq \text{measured body temperature} \leq 37.9^{\circ}\text{C}$ : low fever measured  
body temperature  $\geq 38^{\circ}\text{C}$ : high fever

Figure 12

Figure 13

Figure 14


Figure 15

Prompt

1. Keep the probe and inner cavity clean before measuring
2. Use the thermometer in a stable temperature environment. If the environment temperature changes much, for example from outdoor to indoor, please wait for about 10 minutes before measuring.

3. Don't measure body temperature after measuring extremely high or low temperature. Please wait for 10minutes before measuring body temperature.
4. When the measured object from the environments that is very different from the test environment, please wait for at least half hour before measuring.
5. Trying to avoid using the thermometer when drying the hair, getting wet in the rain, sweating and putting on cosmetics. Don't measure the body temperature within 30 minutes after doing sports, washing and having dinner.

## 2 Object temperature measurement

- Press the measure button to boot and start the measurement
- Make sure the temperature mode on the screen is the object temperature mode.
- Aim the thermometer probe vertically to the object at the distance, 1- 5 cm . Press measure button, after beep in 0.5s the measurement result shows on the screen(Figure 16).
- The thermometer will power off automatically after finishing the measurement in 30s without operation.

Figure 15

### Prompt

1. The function is measuring objects' surface temperature that doesn't represent the internal temperature of the objects
2. Different materials of emission rate might be different. Please refer to the material emissivity, this product default emissivity is 0.95. Different emissivity can cause the deviation between the

measured temperature and real temperature. Take the stainless steel material for example, the measured temperature will be much lower than the actual temperature, please avoid scald.

### 3 Exceed the temperature measuring range.

The Body temperature mode

When the measured value is less than  $32^{\circ}\text{C}$ , the screen shows Lo (Figure 18) with beep indication four times.


When the measured value is higher than  $42.9^{\circ}\text{C}$ , the screen shows Hi (Figure 19) with beep indication four times.

The object temperature mode:

When the measured value is less than  $0^{\circ}\text{C}$ , the screen shows Lo (Figure 20) with beep indication four times.

When the measured value is higher than  $100^{\circ}\text{C}$ , the screen shows Hi (Figure 21) with beep indication four times.

When the environment temperature is less than  $15^{\circ}\text{C}$  or higher than  $35^{\circ}\text{C}$ , the screen shows Er1 (Figure 22), which means it is not allowed to be measured or the accuracy of measurement result cannot be confirmed.


## REPLACEMENT OF THE BATTERY

Open the battery cover, remove the old battery.

Put in 2 new AAA batteries, beep twice that means batteries are installed Pay attention to the direction of the electrodes (Figure 23) .

## Prompt

1. when not used for a long time, please take out the battery, prevent leakage. It is forbidden to put waste battery in the fire.
2. Please properly cope with the waste batteries and avoid the pollution

## MAINTENANCE AND ATTENTIONS

- Since the device would be used frequently, please keep it clean after usage.
- Please keep the sensors and probes inner cavity clean, otherwise it will affect the measuring accuracy.

### Cleaning methods:

1. The surface cleaning: use a clean soft cloth or cotton swab stick with the medical alcohol or water to wipe the dirt.
2. Sensors and sensor inner cavity cleaning: use clean soft cloth or cotton swab dipped a little medical alcohol gently to wipe the probe at the top of the inner cavity or sensors. Start to wipe after alcohol completely evaporate.
  - Before use, please read the instructions, please make sure the battery has been installed
  - Prohibit the thermometer into any liquid, don't place it for long time in very high or very low temperature environment.
  - Avoid collision, falling and mixing with sharp objects, forbidden to disassemble.

- Do not use in strong electromagnetic interference environment.
- The thermometer is placed the place beyond the child's reach
- Suggest to practice to get familiar with the measurement methods, try not to change the factory settings
- The measurement results can not replace the physician diagnosis.
- No other special maintenance is needed during useful life, when fault happens, please contact the vendor or manufacturer.

## TROUBLESHOOTING

The diagnosis	disposal measures
The screen shows „Lo” or „Hi”	<ol style="list-style-type: none"> <li>1. check the measurement object. Unable to ensure the measurement under the condition of drying hair getting wet with water, sweating and putting on cosmetics.</li> <li>2. check the temperature offset settings</li> <li>3. check the operating environment. The environment change influence on the measurement result. If the environment temperature change is too much, immediately measure low temperature object after finishing measuring ultra high temperature object, test error would occur. In this case, please wait for 20 minutes in the test environment before starting measurement.</li> <li>4. check the measuring distance(1 cm to 5 cm).</li> </ol>
No display or display abnormal	<ol style="list-style-type: none"> <li>1. Check the battery and install again</li> <li>2. Check whether the operation is made in the setting menus.</li> </ol>
Buttons have no response	unload the batteries and load again

No indication sound	check if the indication sound is off
Power off when boot	Check the battery, loading and unloading battery again

	Recycling garbage
	Read the operation manual before use
	Warning
	BF type application device
	Batch number
	Production Date

**Guidance and manufacturer 's declaration-electromagnetic emissions**  
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions GB 4824	Group 1	The device is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions GB 4343.1	Class B	
Harmonic radiation GB 17625.1	N/A	
Voltage fluctuations / flicker emissions GB 17625.2	N/A	

Guidance and manufacturer 's declaration-electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the User of the device should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 2, \pm 4, \pm 6$ kV for Contact discharge $\pm 2, \pm 4, \pm 8$ kV air discharge	$\pm 2, \pm 4, \pm 6$ kV for Contact discharge $\pm 2, \pm 4, \pm 8$ kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	$\pm 2, \text{kV}$ for a.c. power lines $\pm 1, \text{kV}$ for d.c. power lines	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 1, \text{kV}$ line(s) to line(s) $\pm 2, \text{kV}$ line(s) to earth	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations in power supply input lines GB/T 17626.11	<5% UT (>95 dip in UT) for 0.5 cycle	N/A	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery
	40% UT	N/A	
	(60% dip in UT) for 5 cycles	N/A	
	70% UT (30% dip in UT) for 25 cycles	N/A	
	<5% UT (>95% dip in UT) for 5 s	N/A	
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment
NOTE UT is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer 's declaration-electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150kHz - 80MHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the device than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1.2 \sqrt{P}$ $d=1.2 \sqrt{P}$ 80MHz-800MHz $d=2.3 \sqrt{P}$ 800MHz-2.5GHz Here P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. d is the recommended separation distance in meters (m).
RF radiation GB/T 17262.3	3 V/m 80 MHz to 2,5 GHz	N/A	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:
NOTE 1: At 90MHz and 800MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<p>a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.</p> <p>b Over the frequency range 150kHz to 80MHz, field strengths should be less than 3V/m.</p>			

Recommended separation distances between portable and mobile RF communications equipment and the product

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter M		
Rated maximum output power	150kHz to 80MHz $d = 1.2 \sqrt{P}$	80MHz to 800MHz $d = 1.2 \sqrt{P}$	800MHz to 2.5GHz $d = 2.3 \sqrt{P}$
0,01	0,01	12	0,23
0,1	0,1	38	0,73
1	1	2	2,3
10	10	8	7,3
100	100	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated by using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## COMMITMENT TO QUALITY AND AFTER-SALES SERVICE

This product provides 1 year free maintenance service.

Note: The free service does not include the failure and damage due to user's personal reasons or unauthorized disassembling.

Prompt:

Please keep the purchase vouchers to facilitate future maintenance.

## Disposal

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

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